


DELight Collaboration Meeting

Lena Hauswald



Uni Freiburg, Freiburg, 18.06.25



Heat flow in dielectric substrates

- Worst case scenario: LAMCAL works as thermal, not athermal detector
- Finite-Elemente-Methode (FEM) simulation with Elmer to determine time resolution for worst-case scenario

Short and spatially
localized heat pulse
(delta-pulse)





Heat flow in dielectric substrates

- Worst case scenario: LAMCAL works as thermal test rather than thermal detector
- Finite-Elemente-Methode (FEM) simulation resolution for worst-case scenario



Ongoing master thesis:

Experimental validation of simulation results

Short and spatially localized heat pulse (delta-pulse)



Quasiparticle detection



- Typical QP energies ≥ 0.8 meV
- Signal gain: Difference between evaporation and adsorption energy
- Dependent on absorber properties:
 - Silicon: $\sim 10\times$;
 - Sapphire: $\sim 20\times$


Magnetic calorimeters on sapphire substrate

Goal:

- First-time realization of a fully operational MMC on a sapphire substrate


Figures of merit validating full MMC functionality:

- Freezing of persistent current
- Agreement between magnetization measurement and theory expectation




Magnetic calorimeters on sapphire substrate

Measurement setup




Slit to themally decouple



Magnetic calorimeters on sapphire substrate

Results:


- First operating MMCs on sapphire
- Persistent current lower then anticipated
 - Dirty substrate, therefore poor Nb quality
- Rising deviation for higher currents and low temperatures
 - Au:Er concentration much higher than intended (Effect of spin-glass?)




Magnetic calorimeters on sapphire substrate

Next to come:

- Fabrication of MMCs with reduced AuEr-Concentration and on clean substrates
- Optimization of parameters for improved freezing of persistent current
- Pulse measurements



Athermal phonon detector



Creation of
athermal phonons



Cooper pair breaking;
Quasiparticles (QPs)
creation



Diffusion


Thermalization through
electron-electron interactions
within temperature sensor

Quasiparticle lifetime experiment

Experiment motivated by
CRESST et al.

Angloher et al., *J Low Temp Phys* 184, 323-329 (2016)

- Determination of diffusion length and quasiparticle lifetime
- Dependent on the choice of material and the quality of material
- Phonon collector geometry optimization





Quasiparticle lifetime experiment



1D Diffusion equation

$$\frac{\delta n}{\delta t} = D \frac{\delta^2 n}{\delta x^2} - \frac{n}{\tau_{QP}}$$



n : Number of QPs
 D : Diffusion constant
 τ_{QP} : QP lifetime




Quasiparticle lifetime experiment



Quasiparticle lifetime experiment



Quasiparticle lifetime experiment




Quasiparticle lifetime experiment

- First pulse measurements with $\text{Kr}^{83\text{m}}$ -source
- Al foil to shield electrons
- Only one operating MMC
 - First pulse shape analysis


Preliminary results following

Quasiparticle lifetime experiment

Slow pulse




Fast pulse




Quasiparticle lifetime experiment

- Rise time estimation (10 – 90 % of amplitude):

Slow pulse




Fast pulse




Quasiparticle lifetime experiment

- Rise time estimation (10 – 90 % of amplitude):

Slow pulse




Fast pulse





Quasiparticle lifetime experiment

- Decay time estimation through exponential fit:



Quasiparticle lifetime experiment

■ Decay time estimation through exponential fit:



Quasiparticle lifetime experiment

■ Decay time estimation through exponential fit:

